

PRODUCT	HIGH VACUUM							PORTABLE	
	HI VAC PLUS 90 Ltr.	HI VAC PLUS 60 Ltr.	HI VAC MS/SS CABINET	HI VAC FOR LIPO SURGERY	HI VAC JR.	EUROVAC ON TROLLEY	MB-36	EUROVAC AC	INSTAD AC
Housing	Plastic Moulded on two wheels	Plastic Moulded on two wheels	MS / Stainless Steel on four castors	Stainless Steel on four castors	MS / Stainless Steel on four castors	ABS Plastic Moulded on four castors	MS Powder Coated, SS Top on two wheels	ABS Plastic Moulded	ABS Plastic Moulded
MAX VACUUM	-700 mm Hg	-600 mm Hg	-700 mm Hg	-700 mm Hg	-650 mm Hg	-575 mm Hg	-700 mm Hg	-575 mm Hg	-700 mm Hg
MAX FLOW	90 Ltrs / min	60 Ltrs / min	50 Ltrs / min	70 Ltrs / min	42 Ltrs / min	45 Ltrs / min	25 Ltrs / min	45 Ltrs / min	25 Ltrs / min
STANDARD OPERATION (voltage)	220 V AC 50 / 60 Hz 580 Watt	220 V AC 50 / 60 Hz 180 Watt	220 V AC 50 / 60 Hz 380 Watt	220 V AC 50 / 60 Hz 380 Watt	220 V AC 50 / 60 Hz 130 Watt	220 V AC 50 / 60 Hz 90 Watt	220 V AC 50 / 60 Hz 180 Watt	220 V AC 50 / 60 Hz 90 Watt	220 V AC 50 / 60 Hz 65 Watt
JARS	2	2	2	2	2	2	2	1	1
NET WEIGHT	23 Kg	23 Kg	37 Kg	36 Kg	15 Kg	8 Kg + 9 Kg Trolley	24 Kg	8 Kg	4 Kg

PRODUCT	PORTABLE		BATTERY OPERATED				MANUAL OPERATED		ELECTRICAL / MANUAL
	INSTAD DC	EUROLITE	EUROVAC AC/DC	INSTAD AC/DC	EUROVAC A	MULTIVAC	P-9	HANDHELD	EP-30
Housing	ABS Plastic Moulded	ABS Plastic Moulded	ABS Plastic Moulded	ABS Plastic Moulded	MS Powder Coated	MS Powder Coated	Plastic Moulded	Plastic Moulded	MS Powder Coated, SS Top on two wheels
MAX VACUUM	-600 mm Hg	-600 mm Hg	-600 mm Hg	-600 mm Hg	-600 mm Hg	-600 mm Hg	-600 mm Hg	-450 mm Hg	Electric -600 mm Hg Pedal -600 mm Hg
MAX FLOW	21 Ltrs / min	12 Ltrs / min	55 Ltrs / min	21 Ltrs / min	14 Ltrs / min	14 Ltrs / min	200 ml / stroke	-	Electric 45 Ltrs / min Pedal 200 ml / stroke
STANDARD OPERATION (voltage)	220 V AC 12 V DC 70 Watt (No Internal Battery)	220 V AC 50 / 60 Hz 90 Watt	220 V AC 50 / 60 Hz 70 Watt 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp)	220 V AC, 50 / 60 Hz, 70 Watt, 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp)	220 V AC 50 / 60 Hz 55 Watt 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp)	220 V AC 50 / 60 Hz 55 Watt 12 V DC Internal Rechargeable Battery 12 V 4.5 Amp	-	-	220 V AC 50 / 60 Hz 90 Watt
JARS	1	1	1	1	1	1	1	1	2
NET WEIGHT	4 Kg	4.5 Kg	11 Kg	6 Kg	5.5 Kg	8.5 Kg	2.5 Kg	192 gm	15.5 Kg

### Our Certifications



### Our Certifications

