

INSTAD

Main Applications

- Hospital Crash Carts
- Patient Transfer
- Long Term Home Care
- Tracheostomy
- EMS / Ambulance
- Bed Side

INSTAD DC




- EMS / Ambulance
- Bed Side



INSTAD AC



TECHNICAL FEATURES

Model	INSTAD AC / DC BATTERY OPERATED SUCTION UNIT	INSTAD DC BATTERY OPERATED (VEHICLE) SUCTION UNIT	INSTAD AC
Housing	ABS Plastic Moulded & 3 Stage LED Indicator for Battery Status	ABS Plastic Moulded	ABS Plastic Moulded
Max. Vacuum / LPM	-760 mm Hg, 30 Ltrs / min	-600 mm Hg, 21 Ltrs / min	-700 mm Hg, 25 Ltrs / min
Pump Type	Piston Pump	Piston Pump	Piston Pump
Jars	1 x 1 Ltr. PC Jar	1 x 1 Ltr. PC Jar	1 x 1 Ltr. PC Jar
Filter	Bacterial Filter Autoclavable / Reusable	Bacterial Filter Autoclavable / Reusable	Bacterial Filter Autoclavable / Reusable
Tubing	8 mm ID x 2 mtr. (Silicon)	8 mm ID x 2 mtr. (Silicon)	8 mm ID x 2 mtr. (Silicon)
Vacuum Gauge	5 cm dia 0-760 mm Hg	5 cm dia 0-760 mm Hg	5 cm dia 0-760 mm Hg
Power Supply	220 V AC, 50 / 60 Hz, 70 Watt 12 V DC Internal Rechargeable Battery (12 V 4.5 Amp) with backup of 150 minutes minimum	220 V AC, 70 Watt, 12 V DC No Internal Battery	220 V AC, 50 / 60 Hz, 65 Watt
Battery	Lead Acid (Ni-cd Optional)	N/A	N/A
Noise Level	55 dB A ± 3	55 dB A ± 3	55 dB A ± 3
Dimensions	39 x 22 x 22 cms	39 x 22 x 22 cms	39 x 22 x 22 cms
Net Weight	6 Kg	4 Kg	4 Kg
Overflow Safety	Mechanical	Mechanical	Mechanical
Accessories / Optional	Car Charger & Power Adaptor (Included) 	Car Charger & Power Adaptor (Included) 	1 Ltr Disposable Jar with liner 110 V on request 
Certifications	CE, IEC 60601-1	CE, IEC 60601-1	CE, IEC 60601-1

ANAND MEDICAIDS PRIVATE LIMITED

Regd. Office:

No.16, Road No.33, Punjabi Bagh Ext. New Delhi-110026, INDIA
 Phone : 91-11-25225225, 25229206
 Fax : 91-11-25225062
 Email : sales@anandind.com

Works:

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA
 Phone: 91-1276-267061

www.anandind.com