

## HANDHELD

### Main Applications

- Patient Transfer
- EMS / Ambulance
- Home Care

## P-9

### Main Applications

- Patient Transfer
- Emergency Department
- Home Care



## TECHNICAL FEATURES

Model	P-9 MANUAL SUCTION UNIT	HAND HELD MANUAL SUCTION UNIT
Housing	Plastic Moulded	Plastic Moulded
Max. Vacuum / LPM	600 mm Hg, 200 ml / push stroke	550 mm Hg
Pump Type	Piston Pump	-
Jars	1 x 1 Ltr. PC Jar	1 x 200 ml Jar
Filter	Bacterial Filter Autoclavable / Reusable	-
Tubing	8 mm ID x 2 mtr. (Silicon)	-
Vacuum Gauge	5 cm dia 0-760 mm Hg	-
Power Supply	N/A	N/A
Noise Level	N/A	N/A
Dimensions	33 x 18 x 30 cms	17 x 5 x 8 cms
Net Weight	2.5 Kg	192 gm
Overflow Safety	Mechanical	N/A
Accessories	Tool Kit, Filter Paper & Push-in Nozzle	Rigid Cannulae & Catheter
Certifications	CE, IEC 60601-1-2012	CE, IEC 60601-1-2012

## ANAND MEDICAIDS PRIVATE LIMITED

### Regd. Office:

No.16, Road No.33, Punjabi Bagh Ext. New Delhi-110026, INDIA  
 Phone : 91-11-25225225, 25229206  
 Fax : 91-11-25225062  
 Email : sales@anandind.com

### Works:

1460 M.I.E, Bahadurgarh, Haryana-124507, INDIA  
 Phone: 91-1276-267061